

# **ChIP Protocol**

Wende/Fong/Duncan  
Kelly Lab  
Revised May 29, 2007

## **DAY ONE (for cells):**

### **Crosslink and cell harvest (C<sub>2</sub>C<sub>12</sub> myotubes):**

1. 2X10 cm plates per condition.
2. Wash plates with PBS (2x10ml).
3. Refresh media 10 ml per plate.
4. Crosslink by adding 108 $\lambda$  37% formaldehyde to each plate.
5. Incubate 10' @ RT on platform shaker.
6. Stop reaction by adding 505 $\lambda$  2.5 M glycine.
7. Incubate 5' @ RT on platform shaker.
8. Rinse w/ 5 ml ice cold PBS (2X).
9. Trypsin/EDTA 1 ml/plate 5' @ 37°C.
10. Add 1 ml PBS; 20% FCS.
  - a. 4ml PBS + 1ml FCS
11. Scrape & transfer to 15 ml conical.
12. Rinse 2 x 1ml/plate ice cold PBS.
13. Pellet 8' 1000rpm (IEC) 4°C. 240xg
14. Wash pellet in 1 ml PBS containing 50  $\mu$ g/ml PMSF.
  - a. 4ml + 12 $\lambda$  100mM PMSF
15. Pellet 8' 1000rpm (IEC) 4°C. 600 x G

### **Cell lysis:**

1. Chill 1 ml Cell Lysis Buffer (CLB) on ice for 5'.
2. Aspirate PBS.
3. Gently resuspend pellet in CLB w/ 1.5x packed cell volume
4. Incubate on ice 10'.
5. Pellet 5' 1500 rpm (IEC) 4°C. 600xg
6. Aspirate CLB.

### **Nuclei lysis:**

1. Chill 4 ml Nuclei Lysis Buffer (NLB) on ice for 3minutes; longer may precipitate SDS.
2. Vigorously resuspend pellet in 1 ml NLB.
3. Incubate 2 hr @ 4°C.

**Sonication:**

1. Chill 13 ml IP Dilution Buffer (IDP) on ice for 5'.
2. Transfer lysed nuclei to 15 ml Falcon 2096 conical tube.
3. Add 800 $\mu$ l IDP.
4. Sonicate on ice.
  - a. 8 repeats
  - b. 30 bursts each
  - c. Output 10%
  - d. Duty 30%
  - e. 1 minute rest between repeats
5. Transfer to microcentrifuge tubes (2ml).
6. Clear insoluble material 16,000 rpm (IEC) 10min @ 4°C microcentrifuge 16Kxg
7. Transfer supernatant to 15ml conical.
8. Dilute to required volume keeping NLB:IDP @ 1:4.
  - a. Vol. = 900 $\mu$ l \* Condition + 300 $\mu$ l Input
  - b. Add ~3.0ml IDP

**DAY ONE (for minced tissue):****Solutions and Buffers to prepare the day of the experiment:**

bracketed volumes are for 2 groups of 3 hearts

double-check volumes needed per experiment!

1. Prepare the following at RT:
  - Butyric Acid (0.5M) [1 mL]
  - PMSF (0.1M) [1 mL]
  - 25X Complete (25X) [4 mL, two tablets]
2. Make the appropriate additions of above solutions to:
  - Cell Lysis Buffer [10mL]
  - Nuclear Lysis Buffer [5mL]
  - IP Dilution Buffer [10mL]

Keep at RT until protocol directs to put on ice.
3. Prepare PBS + Protease Inhibitors:
  - 50 mL 1X PBS
  - 150  $\mu$ L 0.1M PMSF
  - 2.0 mL 25X Complete

**Tissue Harvest and Crosslinking:**

1. Harvest hearts, trimming off atria. Use one adult heart per Ab condition.
2. In a 10cm cell culture dish, mince tissue on ice.
3. Add 10mL of 1X DMEM (4.5g/L glucose) to the dish; keep on ice while other hearts are harvested and minced.
4. Crosslink by adding 270uL of 37% formaldehyde [in the hood!] to each plate.
5. Incubate for 10 minutes at RT on orbital shaker, speed ~50.
6. Stop the reaction by adding 505uL of 2.5M glycine.
7. Incubate for 5 minutes at RT on orbital shaker, speed ~50.
8. Wash tissue 2X with 10mL 1X PBS + protease inhibitors.  
-pipet off as much PBS as possible
9. Scrape and transfer tissue to a 15mL conical using a metal spatula.
10. Rinse dish 1X with 1mL 1X PBS + protease inhibitors and transfer to conical, try to rinse all tissue down to the bottom of the conical.
11. Pellet tissue: spin in IEC 1000rpm for 8 minutes at 4°.
12. Wash each pellet in 1mL 1X PBS + protease inhibitors, try to rinse all tissue down to the bottom of the conical.
13. Pellet tissue: spin in IEC 1000rpm for 8 minutes at 4°.

**Cell Lysis:**

1. Chill Cell Lysis Buffer (CLB) on ice for 5 minutes.
2. Aspirate off PBS.
3. Gently re-suspend tissue pellet in 1.5X packed volume [~1.5 mL] of tissue by pipetting, using a 1000uL tip with the end cut off.
4. Incubate on ice for 10 minutes.
5. Pellet: spin in IEC 1500rpm for 5 minutes at 4°.
6. Aspirate off CLB.

**Nuclei Lysis:**

1. Chill Nuclei Lysis Buffer (NLB) on ice for 3 minutes (do not go longer or SDS may ppt).
2. Re-suspend each pellet in 0.75 mL NLB by pipetting, using a 1000uL tip with the end cut off.
3. Incubate at 4° for 2 hours (in cold room with gentle rocking).  
You will still have visible tissue at the end of this incubation.

**Sonication:**

1. Chill IP Dilution Buffer (IDB) on ice for 5 minutes.
2. Add 600uL of IDB to each lysed nuclei sample.
3. Sonicate on ice:  
40% Duty, Output 5, 10 repeats of 15 seconds.  
Avoid foaming!
4. Transfer samples to 2X 1.5 mL microfuge tubes.
5. Clear insoluble material: spin 15,000 rpm for 10 minutes at 4°.
6. Transfer supernatant to 15 mL conical (recombine samples here)
8. Dilute to required volume, keeping NLB:IDB at 1:4.  
To 0.75 mL NLB + 0.6 mL IDB  
Add 2.4 mL IDB
9. Transfer a 50 uL aliquot of each sample to a 1.5 mL tube (for reverse cross-link).
10. Proceed to Pre-clearing and Immunoprecipitation or snap freeze both aliquots in N<sub>2</sub>(l) and store at -80°.

**THE FOLLOWING IS THE SAME FOR CELLS OR TISSUE:****Pre-clear:**

1. Add 50µl preimmune rabbit serum.
2. Mix @ 4°C for 10min.
3. Add 200µl protein A-Sepharose, 50% slurry.
4. Mix 10min @ 4°C.

**Immunoprecipitation:**

1. Pellet mixture.
  - a. 2000 rpm (IEC) microfuge 2min @ 4°C
2. Aliquot supernatant
  - a. 180λ for input (set @ 4°C until reverse cross-link)
  - b. 900λ for each condition
3. Add antibody (varies considerably).
4. Mix o/n @ 4°C. **END OF DAY 1**
5. Clear insoluble material 16,000 rpm (IEC) microcentrifuge 10min @ 4°C 16Kxg
6. Pipette 60µl Protein-A sepharose/ reaction.
7. Transfer samples to sepharose mix >2 hrs @ 4°C

### **Washing and Elution**

1. Pellet 4000 rpm (IEC) microcentrifuge 2min @ RT-remove to ice **IMMEDIATELY!**
  - a. All washes done with 400 $\lambda$  ice-cold buffers.
  - b. Wash 4x IP wash buffer #1 (20ml)
  - c. Wash 1x IP wash buffer #2 (5ml)
  - d. Wash 3x TE buffer (15ml)
2. Pellet 4000xg microcentrifuge 2' @ RT remove to ice **IMMEDIATELY!**
3. Let sit 1' and remove supernatant.

### **Elution**

4. Add 150 $\mu$ l RT elution buffer. (5ml)
5. Vortex #3 for 15'.
6. Pellet 4500xg 2' @ RT.
7. Transfer supernatant to new tube.
8. Repeat above and combine elutions
  - a. Transfer 140 $\lambda$  after each.

### **Reverse crosslinking:**

1. Add to IP sample:
  - a. 1 $\lambda$  20 mg/ml RNase A and 18 $\lambda$  5M NaCl
2. Add to input sample:
  - a. 1 $\lambda$  20 mg/ml RNase A and 11.5 $\lambda$  5M NaCl
3. Incubate @ 67°C  $\geq$  o/n. **END OF DAY 2**
4. Add 3 $\lambda$  20 mg/ml Proteinase K; gently mix.
5. Incubate @ 45°C for 3 hr.

### **Nucleic acid purification:**

1. Add 120 $\lambda$  TE to input sample; 20 $\lambda$  to IP.
2. Add 1 $\lambda$  10mg/ml tRNA.
3. Extract w/ 300 $\lambda$  1:1 Phenol/Chloroform
  - a. Vortex 10-15"
  - b. Spin full speed 5' RT.
4. Repeat Phenol/Chloroform extraction.
5. Extract w/ 300 $\lambda$  Chloroform
6. Final transfer 300 $\lambda$ .
7. Add 30 $\lambda$  3M NaOAc (pH 5.2), 0.5 $\lambda$  10 mg/ml tRNA, 0.5 $\lambda$  10 mg/ml glycogen.
8. Add 770 $\lambda$  100% EtOH (~2.5X), invert to mix.
9. Incubate 15' @ -80°C then -20°C o/n. **END OF DAY 3**
10. Spin full speed 20' @ 4°C.

11. Wash pellet 800 $\lambda$  70% cold EtOH.
12. Spin full speed 20' @ 4°C.
13. Air dry pellet 5-10' @ RT.
14. Dissolve IP in 30 $\lambda$  & input in 66.7 $\lambda$  TE, pH 8.0.
15. Incubate 5' 37°C, vortex, pulse spin.
16. Store samples @ -20°C. **END OF DAY 4**

**NOTE: Please refer to above protocols (variations exist for most of these steps, specific to cell-type, etc.) and may be altered (incubation times, stopping points, etc.).**